# Hand Therapy Kit

## EXERCISE GUIDE

- Warm up for 3-5 minutes before each session.
- Complete 1-3 sets of 8-12 activity patterns of each exercise.
- Rest approximately 30-60 seconds between each exercise.
- Perform each exercise in a controlled manner.
- If unable to complete 8 activity patterns, increase rest time or decrease number of activity patterns. Or if unable to achieve moderate to maximal fatigue following the completion of 12 activity patterns, decrease rest time or increase number of activity patterns.
- Perform each exercise a minimum of 3 times per week for maximum results.
- Allow 24-48 hours of complete rest between each exercise.

BEFORE BEGINNING THIS OR ANY OTHER EXERCISE PROGRAM, YOU SHOULD ALWAYS CONSULT WITH YOUR DOCTOR OR PHYSICIAN.

### BALL GRIP (WRIST DOWN)

Grip ball in palm. Turn hand in so it faces down. Squeeze the ball in the palm. Hold and relax. Repeat.

### BALL GRIP (WRIST UP)

Grip ball in palm. Turn hand out so it faces up. Squeeze the ball with your palm. Hold and relax. Repeat.

### PINCH

Place ball between the thumb and index finger. Squeeze together. Hold and relax. Repeat.

### THUMB EXTEND

Place ball between bent thumb and two fingers of the same hand. Rolling the ball, extend and straighten the thumb. Repeat.

### OPPOSITION


Place ball in the palm of the hand. Keep between thumb and index finger being exercised. Squeeze thumb and index finger together. Hold and relax hand. Repeat.

### EXTEND OUT

Place the ball on a table. Put tips of fingers on the ball. Roll the ball outward on the table. Repeat.

### SIDE SQUEEZE

Lace ball between any two fingers. Squeeze the two fingers together. Hold and relax. Repeat.

### FINGER BEND

Place the ball in the palm of the hand with fingers pressed into the ball. Push fingers into the ball as if you are bending your fingers. Hold and then relax. Repeat.

Consult with your doctor or physical therapist before beginning any self-treatment program. Suitable for ages 13 and over. Keep out of reach of younger children. This product is not intended to cure any medical conditions.

Gaiam and the Flower of Life logo are trademarks of Sequential Holdco, LLC. ©2024 Sequential Holdco, LLC. All rights reserved. Manufactured and distributed under license by Fit for Life LLC, 1000 S McCaslin Blvd. Suite 300, Superior, CO 80027-8693. **MADE IN CHINA** | SKU 06-58276GUIDE

ACTUAL PRODUCT MAY DIFFER SLIGHTLY FROM ONE SHOWN.

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